

## The flexible solution







New flexibility for transanal and transabdominal surgery



### New flexibility for transanal and transabdominal surgery

Flexibility is the key to efficiency. **KeyPort flex** is efficient in every dimension. The flexible silicone tube simplifies introduction, the intelligent combination of reusable and disposable components creates an economic advantage. Compatibility with the Richard Wolf TEM System and the laparoscopic standard instrument set makes the investment future-proof.

The **KeyPort flex** turnkey system is streamlined, cost-effective and offers great performance. The special telescope and the convenient new rinsing sheath ensure clear vision at all times.

#### Flexible, disposable tube

The sophisticated solution for disposable use is cost-effective and provides exceptionally atraumatic access. High-tech silicone adapts to any anatomy – while also offering a high level of stability and robust construction. The ingenious insertion lug helps to ensure fast and precise placement of the tube which is fixed securely in place using a silicone ring. The well-thought-out geometry of the tube is tapered in the middle section and has a silicone ring at the distal end to offer adequate space for precisie surgery and secure positioning in the anal canal and abdominal wall.

#### Versatile sealing insert

The reusable sealing insert permits extremely flexible working with three valve ports for instrument diameters from 5 to 15 mm. An additional luer connector for passive smoke discharge or active smoke evacuation extends the functionality of the sealing insert. An additional connector for permanent, intrarectal pressure measurement provides compatibility with the Richard Wolf TEM System.









### **DuoRotate and InLine**

Freedom of movement meets precision

#### Three-way division – boundless modularity

All components can be combined with each other in our DuoRotate and InLine instrument generation. Various jaw inserts and sheath tubes can be universally combined with the two handle types (**DuoRotate**, **InLine**). The three-part structure permits easy dismantling for optimal reprocessing.

InLine instruments – double curved sheath tube for maximum freedom of movement.

jaw section, sheath tube positionable in

advance.

#### **Double-curve, single InLine**

The sheath tube of the **InLine** instruments has a double-curve. This curve increases the distance to the neighboring instrument — at the same time, the jaw section moves to the center of the endoscopic image. The resulting effect is that your hand and the instrument tip form an imaginary line. The jaw section of the **InLine** instruments can be turned using the adjustment wheel on the handle. The position of the sheath tube in the handle can be selected before the intervention.










The **DuoRotate** instruments are a pioneering innovation. They elevate Single Port Surgery to a new level. The double rotation function always guarantees the right position for reliable exposure, dissection, or retraction. The sheath tube bent at the distal end can be rotated through 360° using a large and easily

operated rotation wheel on the handle. The toggle lever instantly activates the independently controllable jaw rotation and moves the individual jaw section precisely into position — a technical achievement made in Germany! All this is carried out simply using the toggle lever and the rotation wheel: The hand of

the surgeon can remain in a relaxed position at all times. The ergonomic design of our proven ERAGON handle with disengageable ratchet completes the overall concept.



# **DuoRotate and InLine**

### Freedom of movement meets precision


Jaw insert, insul	ated	DuoRotate sheath tube	InLine sheath tube	Har	ndle	Complete
JL*						instrument  Jaw insert  Sheath tube  Handle  Rinsing adapter (8399300)
		ø 5.8 mm WL 400 mm	ø 5.8 mm WL 400 mm	DuoRotate handle	InLine handle	
	Туре	8399310	8399311	8399301	8399302	Туре
DuoRotate / InLine "Mixter"		•		•		83993010
JL = 21 mm	8399410		•		•	83994010
DuoRotate / InLine "Maryland"		•		•		83993012
JL = 19 mm	8399412		•		•	83994012
DuoRotate / InLine "Metzenbaum"	0200412	•		•		83993013
JL = 18 mm	8399413		•		•	83994013
DuoRotate / InLine "Grasping forceps"	8399414	•		•		83993011
JL = 17 mm	0000414		•		•	83994011
DuoRotate / InLine "Intestine grasping forceps"	8399415	•		•		83993015
JL = 23 mm	0033410		•		•	83994015

<sup>\*</sup>JL = Jaw length when jaws closed (approximate measurements)



Electrode insert, insulated		DuoRotate sheath tube	InLine sheath tube	Handle	Complete instrument  Jaw insert  Sheath tube
					• Handle
		ø 5.8 mm WL 400 mm	ø 5.8 mm WL 400 mm	DuoRotate/InLine electrode handle	
	Type	8399310	8399311	8399303	Туре
DuoRotate / InLine Hook electrode 90°	8399416	•		•	83993016
			•	•	83994016
DuoRotate / InLine Hook electrode 45°		•		•	83993017
			•	•	83994017
Rinsing adapter					

Rinsing adapter				
	Туре			
Rinsing adapter				
	8399300	•		





# **KeyPort** System

	KeyPort Flex System
	KeyPort flex mount autoclavable, inner lumen 24 mm, interface for KeyPort silicone sealings 8850244/8850245 and KeyPort flex silicone tube 485061
	KeyPort flex silicone tube, sterile (pack of 10), working length: 55 mm, Ø 33 mm, for connection to KeyPort flex adapter 885060
	KeyPort flex Set comprising: KeyPort flex adapter (885060), KeyPort flex silicone tube (pack of 10, 485061), KeyPort silicone sealing insert 2 luer (pack of 3, 8850244), adapter with twist lock (885015), adapter with twist lock (885010), adapter with twist lock (885005), plastic sleeve (8850007), plastic sleeve (8850007), seal cap (pack of 10, 89.04), seal cap (pack of 10, 8932), seal cap (pack of 10, 89.08), seal cap (pack of 10, 89.02), Two-way stopcock (8850501)
	KeyPort silicone sealing insert, 2 luer  (pack of 3), autoclavable, with three valve ports, with 2 luer-lock connectors for CO <sub>2</sub> insufflation and active or passive smoke-evacuation, interface for KeyPort flex adapter 8850608850244
	KeyPort silicone sealing insert, 2 luer Set comprising: KeyPort silicone sealing insert 2 luer (pack of 3, 8850244), adapter with twist lock (885015), adapter with twist lock (885010), adapter with twist lock (885005), plastic sleeve (8850010), plastic sleeve (8850005), plastic sleeve (8850007), seal cap (pack of 10, 89.04), seal cap (pack of 10, 8932), seal cap (pack of 10, 8930), seal cap (pack of 10, 8931), seal cap (pack of 10, 89.08), seal cap (pack of 10, 89.02), Two-way stopcock (8850501)
	KeyPort flex TEM Set comprising: KeyPort flex adapter (885060), KeyPort flex silicone tube (pack of 10, 485061), KeyPort silicone sealing insert TEM (pack of 3, 8850245), adapter with twist lock (885015), adapter with twist lock (885010), adapter with twist lock (885005), plastic sleeve (8850007), plastic sleeve (8850007), seal cap (pack of 10, 89.04), seal cap (pack of 10, 8932), seal cap (pack of 10, 8930), seal cap (pack of 10, 89.08), seal cap (pack of 10, 89.02)
	KeyPort silicone sealing insert TEM  (pack of 3) autoclavable, with three valve ports, with 2 luer-lock connectors for CO <sub>2</sub> insufflation and active or passive smoke-evacuation, additional connector for pressure measurement in conjunction  with the TEM combination unit 2232.201, interface for KeyPort flex adapter 8850608850245
	KeyPort silicone sealing insert TEM Set comprising: KeyPort silicone sealing insert TEM (pack of 3, 8850245), adapter with twist lock (885015), adapter with twist lock (885010), adapter with twist lock (8850005), plastic sleeve (88500010), plastic sleeve (8850007), seal cap (pack of 10, 89.04), seal cap (pack of 10, 8932), seal cap (pack of 10, 8931), seal cap (pack of 10, 89.08), seal cap (pack of 10, 89.02)
-	Rinsing sheath for PANOVIEW telescope for PANOVIEW telescope 89354425, rinsing sheath WL 475 mm, Ø 7 mm
	PANOVIEW telescope 5.7 mm,  0°, WL 500 mm, distortion-free, full-screen image