

Afia Examination Table

Top Selling Range of Examination Tables

Lojer's modern, easy to use examination tables have been designed in close collaboration with healthcare professionals to best serve their needs. To ensure consistently high quality, our tables are made in Finland using the latest manufacturing methods and technology. Our popular Afia product family has been specially designed for demanding applications, where patients are moved whilst on the examination table. The tables adapt to a variety of different examination and treatment scenarios, and a number of accessories from our wide range can be combined to meet different needs."

New models: easier to move, more sturdy with design twin castors.



**Afia 4050
Gynaecological
Examination Table**

**AFIA 4140 Gynaecological
Examination Table**



Afia 4040 General Examination Table

LOJER® For easy care

■ Healthcare | Senior Care | Medical Wellness

Afia Examination Tables	4040	4040E	4050	4140
Width	65 cm	65 cm	65 cm	65 cm
Total length	193 cm	193 cm	193 cm	193 cm
Weight	106 kg	101 kg	110 kg	110 kg
Height adjustment	49-94 cm	50-95 cm	49-94 cm	49-94 cm
Back section length, 68 cm	X	X	X	X
Back section adjustment, 0..78°	X	X	X	X
Middle section length, cm	-	-	87	-
Middle section adjustment	-	-	45°	-
Seat adjustment, gas spring assisted, 15°	-	-	X	-
Safe working load (SWL), 200 kg	X	X	X	X
Adjustable heel supports	-	-	X	X
Stainless steel bowl	-	-	X	X
Foot pedal for height adjustment	-	•	-	-
Hand controller for height adjustment	X	X	X	X
Trendelenburg -adjustment, 17°	•	•	•	•
Battery power	•	•	•	•
360° height adjustment foot bar	X	-	X	X
One accessory rail one accessory clamp	X	-	X	X
Extra accessory rail (40cm, max. 6 per table)	•	-	•	•
Clamp for side accessories	•	-	•	•
Angle-adjustable radial clamp for accessories	•	-	•	•
Arm support	•	-	•	•
Anaesthesia frame	•	-	•	•
Paper roll holder	X	X	X	X
Paper cutter	•	•	•	•
Neck pillow	•	•	•	•
Knee supports	-	-	•	•
Side rail	•	•	•	•
IV-pole	•	-	•	•
SKAI Fashion upholstery (13 options)	X	X	X	X
SKAI Clinic upholstery (3 options)	•	•	•	•
SKAI Conductive upholstery (1 option)	•	•	•	•
Centrally locking castors Ø125 mm	X	-	X	X
Individually locking castors Ø125 mm	-	X	-	-
Centrally locking castors Ø150 mm (adds 2 cm to the height range)	•	-	•	•
Directional wheel	X	-	X	X
Electric system, 24V	X	X	X	X
CE-marking	X	X	X	X
Warranty 2 years	X	X	X	X
Country of origin	Finland	Finland	Finland	Finland

Standard feature = X, Optional feature = •, Not available = -