

AFIA 4060/4062 Gynaecological Examination Chairs

Flexible and Easy Adjustments

The Afia 4060/4062 chair models are easy to use examination chairs that have been developed specifically for gynaecological and urological procedures.

It features a two-part top, flexible and easy to use adjustment options and a good level of equipment as standard, such as a neck pillow and paper roll holder. The chairs are available in 17 SKAI imitation leather upholstery options.

The 4060/4062 features electrical height adjustment as well as seat section and back section adjustment; adjustments are made using a hand controller. The chair can also be used in the Trendelenburg position.

Model 4062 also features adjustable calf supports, arm supports and foot controller.

“The Professional's Choice – Versatile Afia Examination Chairs Have Been Made to Work!”



AFIA 4060 Gynaecological Examination Chair

AFIA 4062 Gynaecological Examination Chair

LOJER® For easy care

Healthcare | Senior Care | Medical Wellness

Afia Gynaecological examination chair	4060	4062
Width	62 cm	65 cm
Length	140 cm	140 cm
Weight	75 kg	75 kg
Height adjustment, electric	59-88 cm	59-88 cm
Back section length	90 cm	90 cm
Back section adjustment, electric -7..+70°	-7..+70°	-7..+70°
Seat section adjustment, electrical	X	X
Seat angle, 0..15°	0..15°	0..15°
Safe working load (SWL), 150 kg	X	X
Electric system, 24V	X	X
Adjustable heel supports	X	X
Stainless steel bowl	X	X
Foot control unit	-	X
Hand control unit	X	X
Arm supports	-	X
Paper roll holder	•	•
Neck pillow	•	•
Calf supports	-	X
SKAI Fashion upholstery (13 options)	X	X
SKAI Clinica upholstery (3 options)	•	•
SKAI Conductive upholstery (1 option)	•	•
Castors Ø75mm x 2, Ø125 mm x 2	X	X
CE-marking	X	X
Warranty 2 years	X	X
Country of Origin	Finland	Finland

Standard feature = X Optional feature = • Not available = -

Accessories



Arm supports



Calf supports



Foot controller