

Uniquely sharp and now  
even more ergonomic!



**ERAGON** *modu*  
*blades*

Modular "Metzenbaum" Scissors  
with single-use blades


# ERAGON<sup>modu</sup> blades

## Modular "Metzenbaum" Scissors with single-use blades

**The Best of Both Worlds.** Our popular modular "Metzenbaum" Scissors provide an effective combination of single-use and reusable components: These scissors offer the outstanding cutting quality of single-use scissors, but the low proportion of disposable material makes these scissors much more cost-effective than single-use scissors. Insulation extending over the joint mechanism also offers additional safety, e.g. for dissection close to organs with monopolar HF current. The reusable shaft and handle guarantees high stability and long life. They also bring about a major reduction in running costs.



ø 5.5 mm



ERAGON<sup>modu</sup> blades

WL	Scissor attachment	Pull rod	Sheath	Handle	Complete instrument
	"Metzenbaum", HF mono, (pack of 10), disposable article, packed singly, sterile		insulated	without lock, rotatable, with HF connector	Scissor attachment + Pull rod + Sheath + Handle
<b>260 mm</b>	4393.040	8392934	8392935	83930073	<b>83920407</b>
<b>330 mm</b>		8393934	8393935		<b>83930407</b>
<b>450 mm</b>		8394934	8394935		<b>83940407</b>

- **Guaranteed best cutting quality:**  
**100% sharp every time**
- **Easy assembly:**  
 Screw on scissor attachment – snap handle into place – ready for use!
- **Extremely cost effective:**  
 High-quality, reusable components with a high level of stability and long service life – low proportion of disposable material and therefore low reduced impact on the environment
- **More safety:**  
 Extended insulation covering the joint mechanism
- **Improved ergonomic design:**  
 The rotatable **ERAGON** handle with HF connector facilitates precise operation for delicate work and maximum surgeon comfort

Specifications subject to change without notice.

B 781\_X.12\_09.1