

Trust
the Original













**Miniature Nephroscope
15 / 18 Fr. by Lahme**

Minimum sheath size – maximum stability

Miniature Nephroscope

15 / 18 Fr. by Lahme

Percutaneous miniature low-pressure continuous-irrigation nephroscope by Lahme
comprising:

Endoscope and accessories	
	PANOVIEW operating endoscope , viewing direction 12°, outer ø 12 Fr., WL 225 mm, with laterally offset eyepiece, 6 Fr. working channel, automatic valve with seal membrane 8968.421 also: Sterile universal sealing valve, pack of 5 4712348
	Continuous irrigation sheath , round, WL 205 mm, with rotatable irrigation ring, distal tip straight, automatic locking mechanism 15 Fr. 8968.001 18 Fr. 8968.011
	Obturator , hollow for sheath 8968.001 8968.101 for sheath 8968.011 8968.111
	Amplatz sheath , 18 Fr., WL 150 mm 8968.041
Access instruments	
	Dilator , WL 320 mm, for "one-step" dilatation, stainless steel 12 Fr., suitable for 15 Fr. sheath 8968.012 15 Fr., suitable for 18 Fr. sheath 8968.015
Instruments	
	Stone grasping forceps , ø 2.0 mm, WL 365 mm, with overload protection and mouse tooth geometry 8968.601
	Three-pronged stone grasper , ø 2.0 mm, WL 365 mm 8968.671 also: Spare stone grasper, three-pronged 8968.675
	Ultrasound sonotrode , ø 1.9 mm, WL 344 mm 8962.519
	Grasping forceps "HySafe" , with alligator jaws, 6 Fr., WL 410 mm 8968.6802
	Lithoclast probe , ø 2.0 mm, WL 425 mm 8742.004

- First Mini Nephroscope specially developed for Mini-PCNL: Trust the Original
- For large lower calix stones as an alternative to ESWL
- Also suitable for children and teenagers
- Stainless steel dilators for time-saving "one-step" dilatation
- Automatic locking mechanism between sheath and telescope
- Fast and safe introduction of auxiliary instruments due to automatic seal on working channel
- Smallest possible sheath size with robust construction due to the use of titanium and stainless steel

You can also ask about our matching disposable sets for PCNL.