

# The Future of Diagnostic Thoracoscopy



## Mini Thoracoscopy Set

for Single Port Technique

# Mini Thoracoscopy Set

## for Single Port Technique

This all-in-one Instrument Set for Single Puncture Thoracoscopy allows minimally invasive diagnostic treatment with one small incision. With the small trocar the procedure is possible under local anesthetics and without painful stress for the patient after surgery.

- Very small Thoracoscopy Set with only 5.5 mm diameter
- Optimum image quality with 50.000 pixels
- Large working channel for 3.5 mm instruments
- Very atraumatic method for the patient - reduced pain
- Suction of fluids possible with instrument inserted
- No conventional rod lenses means lens breakage minimized
- The biopsy forceps provide high quality tissue samples
- The working instruments are very easily changed and are always in the field of view
- Ergonomic design



### Operating Laparoscope Set

#### "Leipzig Egg",

5.5 mm Operating Laparoscope with 3.5 mm working channel, direction of view 0°, WL 215 mm, silicate image guide (8920.401), seal (8920.311), seal cap (15176.020) .....8920.4011

#### Flexible Trocar Sleeve, ø 5.5 mm,

Plastic sheath with thread, distal straight, WL 60 mm .....8905.331  
also:

**Trocar**, spherical .....8905.3315

#### Flexible Trocar Sleeve, ø 5.5 mm,

Plastic, WL 75 mm .....8906.051  
also:

**Trocar**, Tip tapered blunt .....8906.151

#### Sample Excision Forceps,

modular, ø 3.5 mm, spoon-shaped, WL 310 mm, insulated, comprising: 8391.406, 8393.0003 .....8391.4063

#### Hook Electrode,

monopolar, ø 3.5 mm, WL 310 mm .....8379.452

#### Coagulation Button Electrode,

monopolar, ø 3.5 mm, WL 310 mm .....8379.462

#### Probe Rod,

ø 3.5 mm, graduated, WL 310 mm .....8379.672

#### Suction Tube,

ø 3 mm, graduated, WL 450 mm .....8380.68

#### Powder spray,

comprising: 8294.03, 8294.82, 8294.60, 2148.95, 127.00, 8291.33

#### Fiber light guide,

2.3 m long .....8061.253

#### Reprocessing Basket

for machine reprocessing and sterilization, for operating laparoscope 8920.4011, Dimensions (W x L x H): 131.5 x 471.5 x 74 mm .....38044.211

**Mini Thoracoscopy Set,  
ø 5.5 mm, complete ..... 8920500**