

ScanAfia X-ICU bed for Intensive Care Units

NEW! Versatile ICU-bed with extensive features

The quality ScanAfia X-ICU bed is a safe choice for hospital intensive care units. Practical design features provide excellent quality-price ratio. Scan Afia X-ICU bed also fulfils the requirements of the new EN 60601-2-52:2010 standard for medical beds.

Although it can take loads up to 340 kg, the ScanAfia X-ICU bed is still light and easy to manoeuvre, thanks to its Ø150mm design twin castors. As a unique feature you may also have electrically lockable set of castors.

The design of the bed emphasizes both ergonomic working and patient comfort and accessibility; its exceptionally low minimum height makes work easier for nurses and makes getting out of bed easier for patients. The back section slide mechanism also improves patient comfort. The bed also features an extremely broad electric Trendelenburg/anti-Trendelenburg adjustment range.

ScanAfia X-ICU's standard features also include:

nurse control panel, electric CPR function, back section slide, battery function and extensive range of quality accessories.

Our many years' experience in designing and manufacturing hospital and ICU beds ensures that the ScanAfia X-ICU provides a long service life with minimal need for maintenance.

NEW!
Practical design with extensive height range, adequate SWL.



ScanAfia X-ICU bed with pair of side rails, electrically locking set of castors, Hygiene Mattress Mediflex

LOJER® For easy care

Scan Afia X-ICU models	ScanAfia X-ICU480	ScanAfia X-ICU490
Length, cm	215	215
Width, 90 cm (incl. Siderails)	90 cm	100 cm
Height adjustment (* depending on the castor, pls see matrix below)	30..82	30..82
Weight, kg	120 kg	122 kg
Recommended mattress size, cm	205 x 78/88 x 13	205 x 78/88 x 13
Safe working load (SWL)	340 kg	340 kg
Centrally locking Ø150mm standard castors	X	X
Centrally locking Ø150mm design twin castors	•	•
Centrally locking Ø125mm castors, electric	•	•
Electric system, 24 V	X	X
Back section adjustment, 0° - 70°	X	X
Foot section angle, 0° - 30°	X	X
Knee break angle, 0° - 18°	X	X
Trendelenburg/ Anti-Trendelengburg, 12°/20°, electric	X	X
Bed lenghtening, 15cm	X	X
Lifting pole	•	•
Lift support	•	•
Urine bag holder	•	•
5th castor	•	•
Bedding stand, foldable	•	•
Monitor stand	•	•
Nurse control panel	X	X
Meal tray	•	•
Handset	X	X
Manual emergency lowering of back section	X	X
Electric CPR function	X	X
Manual (quick) CPR function	•	•
Back section slide	X	X
Breathable ABS plastic base with mattress holders	X	X
Bumper rollers	X	X
Battery	X	X
Side rails, chromed, normal or long version	•	•
Side rails, painted, normal or long version	•	•
Bed ends, in different colors	•	•
CE-marked	X	X
Warranty, 2/ 10 years	X	X
Country of manufacture	Finland	Finland
* Height adjustment range with standard 150 mm castors	34-82 cm	34-82 cm
Height adjustment range with design twin 125 mm castors (electric)	30-78 cm (option)	30-78 cm (option)
Height adjustment range with design twin 150 mm castors	32-80 cm (option)	32-80 cm (option)

Standard feature = X Optional feature = • Not available = -