









# Since 1935

## CHIRANA - European quality



-  Plastic parts, electronic CHIRANA, micromotors CHIRANA, metallic parts
-  Amalgam separator, raw material, machine equipment,
-  Micromotors, scalers, upholstery for patient chair, raw material, machine equipment
-  OP lights, curing lamp, valves, syringes, suction systems
-  Driving units for patient chairs
-  Machine equipment
-  Aluminium casting parts
-  Scalers
-  Scalers

**Made in Slovakia**



# CHIRANA CHEESE E



# Dentist tables

Dentist table 6N



Dentist table 4N



## Mechanical or pneumatic brake of pantographic arm



Pneumatic brake



Mechanical brake



# System of instruments tubes leading

Whip system




Hanging system



# Control panels

## Button-control panel



- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control ( by foot control or manual from keyboard )
- Disinfection programm - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Torque of MM adjustment with automatic stop when set torque is exceeded
- Indication for instrument treatment with oil Smiloil
- free function ( bell – open door )

## Touch-screen control panel



- + Choice of gear ratios of micromotor ( only for brushless MM )
- Program selection (only for brushless micromotor)

## Additional tray for dentist table



### **Tray table under dentist table**

- 1x silicone tray
- 1x stainless tray
- 2x silicone tray
- 2x stainless tray



# Spittoon block

## Suction systems:

- Ejector saliva
- Ejector suction
- CATTANI
- DURR MSBV
- DURR CS1
- DURR CAS1

## Hygiene:

- Clean water bottle
- Disinfection system
- Clean water bottle + disinfection system





# Arms for assistant table

Standard assistant arm 330 mm



Arm with adjustable height



Extended arm



# Assistant tables


EASY 2N – two positions



EASY 4N – four positions



## Control panel on assistant table



- Glass filling function
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Switch buttons for public water or bottle

# Optional accessories on column

Tray table



Monitor



Tray table and monitor




# Arms of trays on column

Pantographic arm



Simple arm



# Trays on column

## **Silicone tray**

1x silicone tray  
1x silicone with handle



## **Stainless tray**


1x stainless tray  
1x stainless with handle



2x silicone tray  
2x silicone tray with handle



2x stainless tray  
2x stainless tray with handle



1x or 2x silicone tray with cup holder  
1x or 2x silicone tray with cup holder + flacons holder



1x or 2x stainless tray with cup holder  
1x or 2x stainless tray with cup holders + flacons holder



# OP lights

## **FARO ALYA LED**

adjustable light intensity 3,000 - 50,000 lux



## **FARO MAIA LED**

adjustable light intensity 3,000 - 35,000 lux



## **LED - A**

adjustable light intensity 7,000 - 38,000 lux



# Foot controls

## **Multifunctional foot control FCR1-9** ( OP light and water on/off )



## **Multifunctional foot control FCR1-7**



## **Multifunctional foot pedal**



## **Foot control with a spray**




## **Foot switch** ( only instruments activation )




# Patient chairs


**SK 1-01**



**SK 1-08**



# Patient chair SK1-01 - upholsteries



Extra cushionary – 10,0 cm thickness




Cushionary - 8, 0 cm thickness



Standard- 6,0 cm thickness



# Dimensions



	A	B	C	D	E	F
Chirana	705	830	23°	18°	--	min.1900/max.2310
Faro Alya	855	830	50°	40°	64	min.1945/max.2355
Faro Maia	855	830	50°	40°	68	min.1945/max.2355
Led-A	605	635	86°	32°	--	min.1915/max.2325

