










Since 1935  
**CHIRANA - European quality**



**Made in Slovakia**



	Plastic parts, electronic CHIRANA, micromotors CHIRANA, metallic parts
	Amalgam separator, raw material, machine equipment,
	Micromotors, scalers, upholstery for patient chair, raw material, machine equipment
	OP lights, curing lamp, valves, syringes, suction systems
	Driving units for patient chairs
	Machine equipment
	Aluminium casting parts
	Scalers
	Scalers

## CHEESE UNI - CHIRANA's news designed for unique solutions

**New design** intended for work with lying patient - in position "12". CHEESE UNI can be equipped from a minimal version suitable for orthodontic rooms up to standard version of treatment unit. CHIRANA CHEESE UNI is designed for left and right handed dentists with rotary dentist and assistant table to right or to left depending on dentist's need.



## CHEESE UNI – from minimum up to maximum



### ... with assistant arm

you can place big and small suction hoses on the assistant table and save 2 places on the dentist table



### ... with tray table

1x silicone or 2x silicone tray table behind the chair



### ... with OP light



### ... with spittoon



### ... with spittoon and OP light





# Dentist table

**5 places:** syringe/ turbine/ micromotor/ scaler / big suction/ small suction



## Control panel

- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control ( by foot control or manual from keyboard )
- Disinfection programm - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control



# Optional suction systems




## Suction systems:

- DURR MSBV valve
- Ready to connect to central wet suction
- Ejector suction



## Hygiene:

- Clean water bottle system
- Disinfection system
- Clean water bottle + disinfection system
- DCI bottles
- Easy to clean Durr filter
- Easy to clean a suction filter
- Easy to clean an oil-catcher



## OP lights

### FARO ALYA LED

adjustable light intensity 3,000 - 50,000 lux



### FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux



### LED - A

adjustable light intensity 7,000 - 38,000 lux



## Foot controls

Multifunctional foot control FCR1-7- WL  
( wire-less )



Multifunctional foot control FCR1-7



Multifunctional foot pedal



Foot control with spray



Foot switch ( only istruments  
activation )



# Patient chair

## SK 1-01

Maximum load of the chair: 150 kg patient + 50 kg additional equipment

Maximum position of the chair:  $870 \pm 20$  mm


Minimum position of the chair:  $420 \pm 20$  mm

Range of the headrest stroke: 165 mm





# Dimensions



	A	B	C	D	E	F
Faro Alya	855	830	50°	40°	64	min.1810/max.2220
Faro Maia	855	830	50°	40°	68	min.1810/max.2220
Led-A	605	635	39°	32°	--	min.1780/max.2190

